ic0934

FOLDABLE EYE
AND TEMPLE MASSAGER

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MASSAGE ACCESSORY
ACCESSOIRE DE MASSAGE

Thank you for purchasing this product. Please read this instruction manual carefully before operating the ic0934 foldable Eye Massager.

Please pay attention to the important safety warnings and keep this manual for future use.

Note: EBI reserves the right to revise the design and description of this product without notice.

CONTENTS

1. Safety warnings	1
2. Environment	2
3. Restrictions	2
4. Maintenance	3
5. Troubleshooting	3
6. Features and benefits	4
7. Components	5
8. Operation	6
8.1. Battery installation	7
8.2. Use	8
8.3. Preset mode	9
8.4. Air pressure switch	9
8.5. Heating switch	10
9. Product specifications	11
10. Limited one (1) year warranty	11
FRANÇAIS	13

1. SAFETY WARNINGS

- This product is not recommended for children.
- Disconnect the massage unit before cleaning.
- Please operate this product according to the instructions In this manual.
- This product is designed for indoor use only.
- Do not exert excessive pressure on the massager. This may significantly slow down or stop the motor, potentially cause overheating and subsequent damage to the massager.
- Press the “stop” button to terminate the massager if discomfort is experienced.
- The total recommended usage time of the massager is 15 minutes.
- Please avoid use if the massage unit is damaged and call your iComfort customer service at: 1-866-967-7333.
- Do not spill liquids onto this product to avoid electric shock or cause damage to this product.

2. ENVIRONMENT

The iComfort massager should be used only in well ventilated, normally dry and dust free environments.

3. RESTRICTIONS

- Remove glasses or contact lens before use.
- Abstain from using the massager in the case of glaucoma, retinoschisis or any serious eye disorders and consult your doctor before using this massage unit again.
- See a physician's advice if you use eyeglasses with a strength of 8.0 or higher.
- People with eye injuries or any other eye health related issues should consult their doctor before using this massage unit.
- Stop immediately the use of the massager if you feel unwell and consult your doctor.

4. MAINTENANCE

- This product should only be maintained by iComfort certified technicians; otherwise, the warranty will be void.
- Make sure the wall outlet is in good working order before use.
- If the massager will not be used for an extended period of time, store in a dry and dustless environment.
- Avoid product exposure to high temperatures and direct sunlight.
- Clean this product with a dry cloth. Do not use thinner, benzene or alcohol.
- We recommend to use the massager intermittently and not for long time.
- Do not insert the batteries if the power adapter is used.
- If the device won't be used for a long time, take out the batteries to avoid battery liquid leakage and causing fault.

5. TROUBLESHOOTING

- It is normal for the motor to make noise while in use.
- If the massager doesn't work normally, please check that the plug and socket are properly connected and the switch is turned on.
- After several use, the massager might power-off because of overheat. It can be used again after half an hour of cooling time.


6. FEATURES AND BENEFITS

Good things really do come in small packages. Fabulously fun and pleasantly powerful, the ic0934 iComfort® eye and temple massager uses air bags to provide a soothing feeling of comfort and relaxation. For extra relief, we've even added a heat feature. The unique foldable ergonomic design help for a perfect fit. The ic0934 eye massager has a 5 minute automatic massage mode. It is also possible to go in manual mode and have the eye massager do one specific massage action at a time. Easy to use, you can appreciate this amazing eye and temple massager anywhere you desire.

- 3 Types of Air Pressure Massage: Ideal for those who simply wish to be pampered in the comfort of their home!
- 2 Types of Heat Massage


7. COMPONENTS


8. OPERATION

Master Unit

Manual Controller

1. Air pressure mode/heating display area
2. Power on/heating selection/ power off key
3. Air pressure switch key

8. 1. BATTERY INSTALLATION

1. Open the battery cover of the manual controller and insert three AA batteries.
2. Close the battery cover

Note: Check whether the polarity + and - of the battery are correctly installed in the unit.

Also check whether the battery cover is closed properly.

8.2. USE

1. Place the Eye Massager over the eyes.

2. Adjust the strap length and place it around the head.

3. Press the ON/OFF button, the device starts to work.


8. 3. PRESET MODE

The default modes are air pressure and heating after power on. The default air pressure mode is Relaxing (1) and the default heating mode is low temperature (h). The display will then display “1”. The default time is 15 minutes.

8. 4. AIR PRESSURE SWITCH

Press the Air pressure switch once after power on and set the air pressure mode to “Kneading”; the air pressure mode/heating display area displays “2”. Press the Air pressure switch twice after power on and set the air pressure mode to “Comfortable”; the air pressure mode/heating display area displays “3”.

Air pressure mode selection chart


8. 5. HEATING SWITCH

Press the 'Power on/heating selection/power off' again after having powered on the unit to set the heating mode to "High temperature". The air pressure mode/heating displays "H" for three seconds indicating the mode is activated.

Press the 'Power on/heating selection/power off' key again to set the heating mode to "Low temperature". The air pressure mode/heating displays "h" for three seconds indicating the mode is activated.

Press the 'Power on/heating selection/power off' key again to disable the heating.

Press and hold the Power on/heating selection/power off key for two seconds to turn off the device.

Heating mode selection chart

9. PRODUCT SPECIFICATIONS

DC:	4.5V 1000mA
Power consumption:	4.5W
Dimensions eye massager:	220mm*100mm* 86mm (L*W*H)
Dimensions controller:	184mm* 36mm* 39mm (L*W*H)
Net weight:	390g

10. LIMITED ONE (1) YEAR WARRANTY

This product is guaranteed against all manufactured faults and defects for a period of one (1) year from the date of purchase. For repairs or product defects, **please call our customer service department at: 1-866-967-7333.** Please have your proof of purchase on hand. If we are not able to replace your defective part, we will send you a service technician for repairs. Important notes: Any damages caused by an accident or an abusive usage are NOT covered under warranty. We also decline all responsibility concerning indirect or consecutive damages. In certain Canadian provinces, it is not permitted to limit indirect or consecutive damages. Therefore it is possible these limitations are not applicable to you. This warranty gives you specific legal rights. You may also have other legal rights which vary from one province to another.

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Pictures are for reference only, colors and models may vary without notice. Made in China. iComfort™ is a trademark of EBI Inc. Photographies pour référence seulement. Les couleurs et modèles sont sujets à varier sans préavis. Fabriqué en Chine. iComfortMC est une marque de commerce de EBI inc.